


Product Name: Recombinant Human FCGR2A Protein
Catalog #: PHH5688



Summary

| | |
|------------------------------|---|
| Name | FCG2A_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P12318 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 20 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | NVLQEDSVTL TCQGARSPE SSIQWFHNGN LIPHTHTQPSY RFKANNND SG EYTCQTGQTS LSDPVHLTVL SEWLVLQTPH LEFQEGETIM LRCHSWKDKP LVKVTFQNG KSQKFSLDLP TFSIPQANHS HSGDYHCTGN IGYTLFSSKP VTITVQVPSM |

image



Background

Alternative Names Low affinity immunoglobulin gamma Fc region receptor II-a (IgG Fc receptor II-a) (CDw32) (Fc-gamma RII-a) (Fc-gamma-RIIa) (FcRII-a) (CD antigen CD32)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FCGR2A Protein
Catalog #: PHH5688

