


Product Name: Recombinant Human HTRA3 Protein
Catalog #: PHH5708



Summary

| | |
|------------------------------|--|
| Name | HTRA3_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P83110 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 15 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | AIPSDRITRF LTEFQDKQIK DWKKRFIGIR MRTITPSLVD ELKASNPDFP EVSSGIYVQE VAPNSPSQRG GIQDGDIVK VNGRPLVDSS ELQEAVLTES PLLLEVRRGN DDLLFSIAPE |

image



Background

Alternative Names Serine protease HTRA3 (EC 3.4.21.-) (High-temperature requirement factor A3) (Pregnancy-related serine protease)

Note

For Research Use Only , Not for Diagnostic Use