



Product Name: Recombinant Human C6orf89 Protein
Catalog #: PHH5711



Summary

Name	CF089_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6UWU4
Host	HEK293
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GVPLHGGDED RFPDFDPWW TNDCEQNESE PIPANCTGCA QKHLKVMLLE DAPRKFERLH PLVIKTGKPL

image



Background

Alternative Names Bombesin receptor-activated protein C6orf89 (Amfion)

Note

For Research Use Only , Not for Diagnostic Use