



Product Name: Recombinant Human RNF130 Protein
Catalog #: PHH5725



Summary

Name	GOLI_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q86XS8
Host	HEK293
Species	Human
Molecular Mass	61 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQEYYTALIN VTVQEPGRGA PLTFRIDRGR YGLDSPKAEV RGQVLAPLPL HGVADHLGCD PQTRFFVPPN IKQWIALLR GNCTFKEKIS RAAFHNAVAV VIYNNKSKEE PVTMTHPGTG DIIAVMITEL RGKDILSYLE KNISVQMTIA VGTRMPPKNF

image



Background

Alternative Names	E3 ubiquitin-protein ligase RNF130 (EC 2.3.2.27) (Goliath homolog) (H-Goliath) (RING finger protein 130) (RING-type E3 ubiquitin transferase RNF130)
--------------------------	--

Note

Product Name: Recombinant Human RNF130 Protein
Catalog #: PHH5725



For Research Use Only , Not for Diagnostic Use