


Product Name: Recombinant Human ERFE Protein
Catalog #: PHH5732




Summary

Name	ERFE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q4G0M1
Host	HEK293
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPGAPSRRA RREPPPGNEL PRGPGESRAG PAARPPEPTA ERAHSVDPRD AWMLFVRQSD KGVNGKKRSR GKAKKLKFLG PGPPGPPGPQ GPPGPIIPPE ALLKEFQLLL KGAVRQRERA

image



SDS



WB

Background

Alternative Names	Erythroferrone (Complement C1q tumor necrosis factor-related protein 15) (Myonectin)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use