



Summary

Name	I17RD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NFM7
Host	HEK293
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQLAVAAGGS GRARGADTCG WRGVGPASRN SGLYNITFKY DNCTTYLNPV GKHVIADAQN ITISQYACHD QVAVTILWSP GALGIEFLKG FRVILEELKS EGRQCQQLIL KDPKQLNSSF KRTGMESQPF LNMKFETDYF VKVVPFPSIK NESNYHPFFF RTRACDLLLQ PDNLACKPFW KPRNLNISQH GSDMQVSFDH APHNFGFRFF YLHYKLNKHEG PFKRKTCKQE QTTETTSCLL QNVSPGDYII ELVDDTNTTR KVMHYALKPV

image



Background

Alternative Names	Interleukin-17 receptor D (IL-17 receptor D) (IL-17RD) (IL17Rhom) (Interleukin-17 receptor-like protein) (Sef homolog) (hSef)
--------------------------	---

Product Name: Recombinant Human IL17RD Protein
Catalog #: PHH5739



Note

For Research Use Only , Not for Diagnostic Use