


Product Name: Recombinant Human BOLA3 Protein
Catalog #: PHH5772




Summary

Name	BOLA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST & C-His
Accession #	Q53S33
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAAWSPAAAA PLLRGIRGLP LHHRMFATQT EGELRVQIL KEKFPRATAI KVTDISGGCG AMYEIKIESE EFKEKRTVQQ HQMVNQALKE EIKEMHGLRI FTSVPKR

image



SDS



WB

Background

Alternative Names BOLA-like protein 3

Note

For Research Use Only , Not for Diagnostic Use