


Product Name: Recombinant Human MCOLN2 Protein
Catalog #: PHH5778



Summary

Name	MCLN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	Q8IZK6
Host	HEK293
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VAFKEDNTVA FKHLFLKGYS GTDEDDYSCS VYTQEDAYES IFFAINQYHQ LKDITLGLTG YGENEDNRIG LKVCKQHYKK GTMFPSNETL NIDNDVELDC VQLDLQDLSK KPPDWKNSSF FRLEFYRLQ VEISFHLKGI DLQTIHSREL PDCYVFQNTI IFDNKAHSGK IKIYFSDSAK IEECKDLNIF

image



Background

Alternative Names Mucolipin-2 (Transient receptor potential channel mucolipin 2) (TRPML2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MCOLN2 Protein
Catalog #: PHH5778

