



Product Name: Recombinant Human CGB3 Protein
Catalog #: PHH5786



Summary

Name	CGB3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P0DN86
Host	HEK293
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SKEPLRPRCR PINATLAVEK EGCPVCITVN TTICAGYCPT MTRVLQGVLP ALPQVVCNYR DVRFESIRLP GCPRGVNPVV SYAVALSCQC ALCRRSTTDC GGPKDHPLTC DDPFRQDSSS SKAPPSLPS PSRLPGPSDT PILPQ

image



Background

Alternative Names Choriogonadotropin subunit beta 3 (Choriogonadotropin subunit beta) (CG-beta) (Chorionic gonadotropin chain beta)

Note

For Research Use Only , Not for Diagnostic Use