



Product Name: Recombinant Human CNTNAP5 Protein
Catalog #: PHH5792



Summary

Name	CNTP5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WYK1
Host	HEK293
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DDPLASLLSP MAFSSSSDLT GTHSPAQLNW RVGTGGWSPA DSNAQQWLQM DLGMRVEITA VATQGRYGSS DWVTSYSLMF SDTGRNWKQY KQEDSIWTFA GNMNADSVVH HKLLHSVRAR FVRFVPLEWN PSGKIGMRVE VYGCSYKSDV ADFDGRSSLL YRFNQKLMST

image



Background

Alternative Names Contactin-associated protein-like 5 (Cell recognition molecule Caspr5)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CNTNAP5 Protein
Catalog #: PHH5792

