


Product Name: Recombinant Human CRISP2 Protein
Catalog #: PHH5794




Summary

Name	CRIS2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P16562
Host	HEK293
Species	Human
Molecular Mass	67 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDPAFTALLT TQLQVQREIV NKHNELRKAV SPPASNMLKM EWSREVTNA QRWANKCTLQ HSDPEDRKTS TRCGENLYMS SDPTSWSSAI QSWYDEILDF VYGVGPKSPN AVVGHYTQLV WYSTYQVGCG IAYCPNQDSL KYYYVCQYCP AGNNMNRKNT PYQQGTPCAG CPDDCDKGLC TNSCQYQDLL SNCDCLKNTA GCEHELLKEK CKATCLCENK IY

image



SDS



WB

Background

Alternative Names	Cysteine-rich secretory protein 2 (CRISP-2) (Cancer/testis antigen 36) (CT36) (Testis-specific protein TPX-1)
--------------------------	--

Note

Product Name: Recombinant Human CRISP2 Protein
Catalog #: PHH5794



For Research Use Only , Not for Diagnostic Use