


Product Name: Recombinant Human CCL27 Protein
Catalog #: PHH5820



Summary

| | |
|------------------------------|--|
| Name | CCL27_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-His & N-SUMO |
| Accession # | Q9Y4X3 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 22 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | FLLPPSTACC TQLYRKPLSD KLLRKVIQVE LQEADGDCHL QAFVLHLAQR SICIHQPQNPS LSQWFEHQER KLHGTLPKLN FGMLRKMG |

image



Background

| | |
|--------------------------|--|
| Alternative Names | C-C motif chemokine 27 (CC chemokine ILC) (Cutaneous T-cell-attracting chemokine) (CTACK) (ESkine) (IL-11 R-alpha-locus chemokine) (Skinkine) (Small-inducible cytokine A27) |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use