



Product Name: Recombinant Human CLCN3 Protein
Catalog #: PHH5823



Summary

Name	CLCN3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P51790
Host	HEK293
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YYRNSYNSIT SASSDEELLD GAGVIMDFQT SEDDNLLDGD TAVGTHYMT NGGSINSSTH LLDLLDEPIP GVGTYDDFHT IDWVREKCKD RERHRRINSK KKESAWEMTK

image



Background

Alternative Names	H(+)/Cl(-) exchange transporter 3 (Chloride channel protein 3) (ClC-3) (Chloride transporter ClC-3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use