


Product Name: Recombinant Human NCR2 Protein
Catalog #: PHH5827



Summary

Name	NCTR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	O95944
Host	HEK293
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VAGQTLTVRC QYPPTGSLYE KKGWCKEASA LVCIRLVTSS KPRTMAWTSR FTIWDDPDAG FFTVTMTDLR EEDSGHYWCR IYRPSDNSVS KSVRFYLVVS PASASTQTSW TPRDLVSSQT QTQSCVPPTA GARQAPESPS TIPVPSQPQN STLRPGPAAP

image



Background

Alternative Names	Natural cytotoxicity triggering receptor 2 (Lymphocyte antigen 95 homolog) (NK cell-activating receptor) (Natural killer cell p44-related protein) (NK-p44) (NKp44) (CD antigen CD336)
--------------------------	--

Note

Product Name: Recombinant Human NCR2 Protein
Catalog #: PHH5827



For Research Use Only , Not for Diagnostic Use