


Product Name: Recombinant Human TMEM30B Protein
Catalog #: PHH5830



Summary

Name	CC50B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q3MIR4
Host	HEK293
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KELEYDYTGD PGTGNCVCA AAGQGRALPP PCSCAWYFSL PELFQGPVYL YYELTNFYQN NRRYGVS RDD AQLSGLPSAL RHPVNECAPY QRSAAGLPIA PCGAIANSLF NDSFSLWHQR QPGGPYVEVP LDRSGIAWWT DYHVKFRNPP LVNGSLALAF QGTAPPPNWR RPVYELSPDP NNTGFINQDF VVWMRTAALP TFRKLYARIR QGNYSAGLPR GAYRVNITYN YPVRAFGG HK LLISSISWM

image



Background

Alternative Names Cell cycle control protein 50B (P4-ATPase flippase complex beta subunit TMEM30B) (Transmembrane protein 30B)

Note

Product Name: Recombinant Human TMEM30B Protein
Catalog #: PHH5830



For Research Use Only , Not for Diagnostic Use