

Product Name: Recombinant Human SLC2A10 Protein
Catalog #: PHH5841



Summary

Name	GTR10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95528
Host	HEK293
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VPMDSGPSCL AVPNATGQTG LPGDSGLLQD SSLPPIPRTN EDQREPILST AKKTKPHPRS GDPSAPPRLA LSSALPGPPL PARGHALLRW

image



Background

Alternative Names Solute carrier family 2, facilitated glucose transporter member 10 (Glucose transporter type 10) (GLUT-10)

Note

For Research Use Only , Not for Diagnostic Use