



Product Name: Recombinant Human ITGA8 Protein
Catalog #: PHH5845



Summary

Name	ITA8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P53708
Host	HEK293
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TEFRDKLSPI NISLNYSLDE STFKEGLEVK PILNYYRENI VSEQAHLVD CGEDNLCVPD LKLSARPDKH QVIIGDENHL MLIINARNEG

image



Background

Alternative Names	Integrin alpha-8 [Cleaved into: Integrin alpha-8 heavy chain; Integrin alpha-8 light chain]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use