


Product Name: Recombinant Human CD47 Protein
Catalog #: PHH5857



Summary

Name	CD47_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q08722
Host	HEK293
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LFNKTKSVEF TFCNDTVVIP CFVTNMEAQN TTEVYVKWKF KGRDIYTFDG ALNKSTVPTD FSSAKIEVSQ LLKGDASLKM DKSDAVSHTG NYTCEVTELT REGETIIELK

image



Background

Alternative Names Leukocyte surface antigen CD47 (Antigenic surface determinant protein OA3) (Integrin-associated protein) (IAP) (Protein MER6) (CD antigen CD47)

Note

For Research Use Only , Not for Diagnostic Use