



Product Name: Recombinant Human BCAP29 Protein
Catalog #: PHH5862



Summary

Name	BAP29_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UHQ4
Host	HEK293
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LSNKGVLKTQ AENTNKAACK FMEENEKLR ILKSHGKDEE CVLEAENKKL VEDQEKLKTE LRKTSDALSK AQNDVMEMKM QSERLSKEYD QLLKEHSELQ DRLERGNKKR

image



Background

Alternative Names B-cell receptor-associated protein 29 (BCR-associated protein 29) (Bap29)

Note

For Research Use Only , Not for Diagnostic Use