


Product Name: Recombinant Human CA10 Protein
Catalog #: PHH5863



Summary

Name	CAH10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NS85
Host	HEK293
Species	Human
Molecular Mass	33 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GWWAYKEVVQ GSFVPVPSFW GLVNSAWNLC SVGKRQSPVN IETSHMIFDP FLTPLRINTG GRKVSGMTMYN TGRHVSLRLD KEHLVNISGG PMTYSHRLEE IRLHFGSEDS QGSEHLLNGQ AFSGEVQLIH YNHELYTNVT EAAKSPNGLV VVSIFIKVSD SSNPFLNRML NRDTITRITY KNDAYLLQGL NIEELYPETS SFITYDGSMIT IPPCYETASW IIMNKPVIIT RMQMHSRLRL SQNQPSQIFL SMSDNFRPVQ PLNNRCIRTN

image



Background

Alternative Names	Carbonic anhydrase-related protein 10 (Carbonic anhydrase-related protein X) (CA-RP X) (CARP X) (Cerebral protein 15)
--------------------------	---

Product Name: Recombinant Human CA10 Protein
Catalog #: PHH5863



Note

For Research Use Only , Not for Diagnostic Use