



Product Name: Recombinant Human FNDC5 Protein
Catalog #: PHH5865



Summary

Name	FNDC5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NAU1
Host	HEK293
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TVRHLKANS A VVSWDVLEDE VVIGFAISQQ KKDVRMLRFI QEVNTTTRSC ALWDLEEDTE YIVHVQAISI QGQSPASEPV LFKTPREA EK MASKNKDEVT MKEMGRNQQL

image



Background

Alternative Names Fibronectin type III domain-containing protein 5 (Fibronectin type III repeat-containing protein 2) [Cleaved into: Irisin]

Note

For Research Use Only , Not for Diagnostic Use