


Product Name: Recombinant Human CTCRC Protein
Catalog #: PHH5870



Summary

| | |
|------------------------------|---|
| Name | CTRC_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q99895 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 29 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | CGVPSFPPNL SARVVGGEDA RPHSWPWQIS LQYLKNDTWR HTC GGTLIAS NFVLTA AHCI SNTRTYRVAV GKNNLEVEDE EGSLFVGVDI IHVHKRWNAL LLRNDIALIK LAEVELSDT IQVACLPEKD SLLPKDYPCY VTGWGRLWTN GPIADKLQQG LQPVVDHATC SRIDWWGFRV KKT MVCAGGD GVISACNGDS GGPLNCQLEN GSWEVFGIVS FGSRRGCNTR KKPVVYTRVS AYIDWINEKM QL |

image



Background

Alternative Names Chymotrypsin-C (EC 3.4.21.2) (Caldecrin)

Note

Product Name: Recombinant Human CTRC Protein
Catalog #: PHH5870



For Research Use Only , Not for Diagnostic Use