



Summary

Name	F187B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q17R55
Host	HEK293
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SISCPGKQC QQALLSGNDI LLYCNSSGAH WYYLFTQGKK GRLTSLTNIS NMEIMPEGSL LIKDPLPSQT GLYHCWNKNG RQVVQYEIDF QDVTTLHITH KDLGQRPLQN ETLHLGSKQL IFTWWEPWQD CNRCEEPGEC

image



Background

Alternative Names Protein FAM187B (Transmembrane protein 162)

Note

For Research Use Only , Not for Diagnostic Use