


Product Name: Recombinant Human FZD4 Protein
Catalog #: PHH5899



Summary

Name	FZD4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9ULV1
Host	HEK293
Species	Human
Molecular Mass	47 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EERRCDPIRI SMCQNLGYNV TKMPNLVGHE LQTDALQLT TFTPLIQYGC SSQLQFFLCS VYVPMCTEKI NIPIGPCGGM CLSVKRRCEP VLKEFGFAWP ESLNCSKFPP QNDHNHMCME GPGDEEVPLP HKTPIQPGEE CHSVGTNSDQ YIWVKRSLNC VLKCGYDAGL

image



Background

Alternative Names Frizzled-4 (Fz-4) (hFz4) (FzE4) (CD antigen CD344)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FZD4 Protein
Catalog #: PHH5899

