


Product Name: Recombinant Human EIF1B Protein
Catalog #: PHH5915




Summary

Name	EIF1B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST & C-His
Accession #	O60739
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	STIQNLQSFD PFADATKGDD LLPAGTEDYI HIRIQQRNGR KTLTTVQGIA DDYDKKKLVK AFKKKFACNG TVIEHPEYGE VIQLQGDQRK NICQFLLEV IVKEEQLKVH GF

image



SDS



WB

Background

Alternative Names	Eukaryotic translation initiation factor 1b (eIF1b) (Protein translation factor SUI1 homolog GC20)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use