



Product Name: Recombinant Human EVI2A Protein
Catalog #: PHH5917



Summary

Name	EV12A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P22794
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TNPITPEVDY KGNSTNMPET SHIVALTSKS EQELYIPSVV SNSPSTVQSI ENTSKSHGEI FKKDVCAENN

image



Background

Alternative Names Protein EVI2A (Ecotropic viral integration site 2A protein homolog) (EVI-2A)

Note

For Research Use Only , Not for Diagnostic Use