


Product Name: Recombinant Human GPR37L1 Protein
Catalog #: PHH5920



Summary

Name	G37L1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	O60883
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QSRSKRGTED EEAKGVQQYV PEEWAEYPRP IHPAGLQPTK PLVATSPNPG KDGGT PDSGQ ELRGNLTGAP

image



Background

Alternative Names G-protein coupled receptor 37-like 1 (Endothelin B receptor-like protein 2) (ETBR-LP-2)

Note

For Research Use Only , Not for Diagnostic Use