



Product Name: Recombinant Human CPO protein
Catalog #: PHH5924



Summary

Name	CBPO_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q8IVL8
Host	HEK293
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NIDGYIYTWT TDRLWRKSRS PHNNGTCFGT DLNRNFNASW CSIGASRNCQ DQTFCGTGPV SEPETKAVAS FIESKKDDIL CFLTMHSYGQ LILTPYGYTK NKSSNHPEMI

image



Background

Alternative Names Carboxypeptidase O (CPO) (EC 3.4.17.-)

Note

For Research Use Only , Not for Diagnostic Use