


Product Name: Recombinant Human NRG3 Protein
Catalog #: PHH5935



Summary

Name	NRG3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P56975
Host	HEK293
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FPKAMETTTT TTSTTSPATP SAGGAASSRT PNRISTRLLT ITRAPTRFPG HRVPIRASPR STTARNTAAP

image



Background

Alternative Names Pro-neuregulin-3, membrane-bound isoform (Pro-NRG3) [Cleaved into: Neuregulin-3 (NRG-3)]

Note

For Research Use Only , Not for Diagnostic Use