


Product Name: Recombinant Human FADS3 Protein
Catalog #: PHH5952




Summary

Name	FADS3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5Q0
Host	HEK293
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGPAQPGAPL PTFCWEQIRA HDQPGDKWLV IERRVYDISR WAQRHPGGSR LIGHHGAEDA TDAFRAFHQD LNFVRKFLQP LLIGELAPEE PSQDGPLNAQ LVEDFRALHQ AAEDMKLFDA

image



SDS



WB

Background

Alternative Names	Fatty acid desaturase 3 (EC 1.14.19.-) (Delta(13) fatty acid desaturase) (Delta(13) desaturase)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use