



Product Name: Recombinant Human CFP Protein
Catalog #: PHH5965



Summary

Name	PROP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P27918
Host	HEK293
Species	Human
Molecular Mass	49 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DPVLCFTQYE ESSGKCKGLL GGGVSVEDCC LNTAFAYQKR SGGLCQPCRS PRWSLWSTWA PCSVTCSEGS QLRYYRCVGVW NGQCSGKVAP GTLEWQLQAC EDQQCCPEMG GWSGWGPWEP CSVTCCKGTR TRRRACNHPA PKCGGHCPGQ AQESEACDTQ QVCPHGAWA TWGPWTPCSA SCHGGPHEPK ETRSRKCSAP EPSQKPPGKP CPGLAYEQRR CTGLPPCPVA GGWGPWGPVS PCPVTCLGQ TMEQRTCNP VPQHGGPFCA GDATRTHICN TAVPCPVDGE WDSWGEWSPC IRRNMKISIC QEIPGQQSRG RTCRGRKFDG HRCAGQQQDI RHCYSIQHCP LKGSWSEWST WGLCMPPCGP NPTRARQRLC TPLLKPYPPT VSMVEGQGEK NVTFWGRPLP RCEELQGQKL VVEEKRPLH VPACKDPEEE EL

image



Background

Product Name: Recombinant Human CFP Protein
Catalog #: PHH5965



Alternative Names

Properdin (Complement factor P)

Note

For Research Use Only , Not for Diagnostic Use