


Product Name: Recombinant Human GIPR Protein
Catalog #: PHH5979



Summary

Name	GIPR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His & N-SUMO
Accession #	P48546
Host	HEK293
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TAGELYQRWE RYRRECQETL AAAEPPSGLA CNGSFDMYVC WDYAAPNATA RASCPWYLPW HHHVAAGFVL RQCGSDGQWG LWRDHTQCEN PEKNEAFLDQ

image



Background

Alternative Names	Gastric inhibitory polypeptide receptor (GIP-R) (Glucose-dependent insulinotropic polypeptide receptor)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use