


Summary

Name	GPC6A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q5T6X5
Host	HEK293
Species	Human
Molecular Mass	61 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ITTIPNVKKI GKVVGFAFRR GNISSFHSFL QNLHLLPSDS HKLLHEYAMH LSACAYVKDT DLSQCIFNHS QRTLAYKANK AIERNFVMRN DFLWDYAEPG LIHSIQLAVF ALGYAIRDLC QARDCQNPNA FQPWELLGVL KNVTFDQGNW SFHFDAHGDL NTGYDVVLWK EINGHMTVTK MAEYDLQNDV FIIPDQETKN EFRNLKQIQS KCSKECSPGQ MKKTTRSQHI CCYECQNCPE NHYTNQTDMP HCLLCNNKTH WAPVRSTMCF EKEVEYLNWN

image



Background

Alternative Names	G-protein coupled receptor family C group 6 member A (hGPCR6A) (G-protein coupled receptor GPCR33) (hGPCR33)
--------------------------	--

Product Name: Recombinant Human GPRC6A Protein
Catalog #: PHH5990



Note

For Research Use Only , Not for Diagnostic Use