



Product Name: Recombinant Human IL18BP Protein
Catalog #: PHH5992



Summary

Name	I18BP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	O95998
Host	HEK293
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TPVSQTTTAA TASVRSTKDP CPSQPPVFPA AKQCPALEVT WPEVEVPLNG TSLSCVACS RFPNFSILYW LGNGSFIEHL PGRLWEGSTS RERGSTGTQL CKALVLEQLT PALHSTNFSC VLVDPEQVVQ RHVVLAQLWA GLRATLPPTQ EALPSSHSSP QQQG

image



Background

Alternative Names Interleukin-18-binding protein (IL-18BP) (Tadekinig-alfa)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human IL18BP Protein
Catalog #: PHH5992

