


Summary

Name	CSF3R_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99062
Host	HEK293
Species	Human
Molecular Mass	68 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VSAPIVHLGD PITASCIKQ NCSHLDPEPQ ILWRLGAELQ PGGRQQRLSD GTQESIITLP HLNHTQAFSL CCLNHWGNSLQ ILDQVELRAG YPPAIPHNL CLMNLTTSSL ICQWEPGPET HLPTSFTLKS FKSARGNCQTQ GDSILDCVPK DGQSHCCIPR KHLLLYQNMG IWVQAENALG TSMSPQLCLD PMDVVKLEPP MLRTMDPSPE AAPPQAGCLQ LCWEPWQPGL HINQKCELRH KPQRGEASWA LVGPLEAL QYELCGLLPA TAYTLQIRCI RWPLPGHWSW WSPSLELRTT ERAPTURLDT WWRQRQLDPR TVQLFWKPVP LEEDSGRIQG YVSWRPSGQ AGAILPLCNT TELSCTFHLP SEAQEVALVA YNSAGTSRPT PVVFESRGP ALTRLHAMAR DPHSLWVGWE PPNPWPQGYV IEWGLGPPSA SNSNKTWRME QNGRATGFLL KENIRPFQLY EIIVTPLYQD TMGPSQHVYA YSQEMAPSHA PELHLKHIGK TWAQLEWVPE PPELGKSPLT HYTIFWTNAQ NQSFSAILNA SSRGFVLHGL EPASLYHIHL MAASQAGATN STVLTLMTLT

image



Product Name: Recombinant Human CSF3R Protein
Catalog #: PHH6001



Background

Alternative Names	Granulocyte colony-stimulating factor receptor (G-CSF receptor) (G-CSF-R) (CD antigen CD114)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use