



Product Name: Recombinant Human LRRC3B Protein
Catalog #: PHH6010



Summary

Name	LRC3B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96PB8
Host	HEK293
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSSGGLNVTC SNANLKEIPR DLPPETVLLY LDSNQITSIP NEIFKDLHQL RVLNLSKNGI EFIDEHAFKG VAETLQTLDL SDNRIQSVHK NAFNNLKARA RIANNPWHCD CTLQQVLRSM ASNHETAHNV ICKTSVLDEH AGRPFLNAAN DADLCNLPPK

image



Background

Alternative Names Leucine-rich repeat-containing protein 3B (Leucine-rich repeat protein 15)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LRRC3B Protein
Catalog #: PHH6010

