


Product Name: Recombinant Human TP53I13 Protein
Catalog #: PHH6018




Summary

Name	P5I13_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q8NBR0
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQGPRQPSSS GAKRRRLRAA LGPQPTRSAL RFPSASPGSL KAKQSMAGIP GRESNAPSVP TVSLPGAPG GNASSRTEAQ VPNGQGSPGG CVCSSQASPA PRAAAPPRAA

image



SDS



WB

Background

Alternative Names	Tumor protein p53-inducible protein 13 (Damage-stimulated cytoplasmic protein 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use