


**Product Name: Recombinant Human PF4V1 Protein**  
**Catalog #: PHH6022**



## Summary

<b>Name</b>	PF4V_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P10720
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	9 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	FARAEAEEDG DLQCLCVKTT SQVRPRHITS LEVIKAGPHC PTAQLIATLK NGRKICLDLQ ALLYKIIKE HLES

## image



## Background

<b>Alternative Names</b>	Platelet factor 4 variant (C-X-C motif chemokine 4 variant) (CXCL4L1) (PF4alt) (PF4var1) [Cleaved into: Platelet factor 4 variant(4-74); Platelet factor 4 variant(5-74); Platelet factor 4 variant(6-74)]
--------------------------	--

## Note

For Research Use Only , Not for Diagnostic Use