



Product Name: Recombinant Human PCDHA11 Protein
Catalog #: PHH6029



Summary

Name	PCDAB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9Y5I1
Host	HEK293
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RFPLEGASDA DIEENALLTY RLSKNEYFSL DSPTNGKQIK RLSLILKKSL DREKPELNL LLTATDGGKP ELTGTVRLLV QVLVDVNDNDP EFDKSEYKVS LMENAAKETL VLKLNATDRD EGVNGEVTYS LMSIKPNGRH LFTLDQNNGE VRVNGTLDYE ENKFYKIEVQ ATDKGTPPMA GHCTVWVEIL

image



Background

Alternative Names Protocadherin alpha-11 (PCDH-alpha-11)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PCDHA11 Protein
Catalog #: PHH6029

