


Product Name: Recombinant Human PRTN3 Protein
Catalog #: PHH6041



Summary

Name	PRTN3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P24158
Host	HEK293
Species	Human
Molecular Mass	66 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IVGGHEAQPH SRPYMASLQM RGNPGSHFCG GTLIHPSFVL TAAHCLRDIP QRLVNVVLGA HNVRTQEPTQ QHFSVAQVFL NNYDAENKLN DVLLIQLSSP ANLSASVATV QLPQQDQVPV HGTQCLAMGW GRVGAHDPPA QVLQELNVTV VTFFCRPHNI CTFVPRRKAG ICFGDSGGPL ICDGIIQGID SFVIWGCATR LFPDFFTRVA LYVDWIRSTL R

image



Background

Alternative Names	Myeloblastin (EC 3.4.21.76) (AGP7) (C-ANCA antigen) (Leukocyte proteinase 3) (PR-3) (PR3) (Neutrophil proteinase 4) (NP-4) (P29) (Wegener autoantigen)
--------------------------	--

Note

Product Name: Recombinant Human PRTN3 Protein
Catalog #: PHH6041



For Research Use Only , Not for Diagnostic Use