


Product Name: Recombinant Human SLC22A10 Protein
Catalog #: PHH6055



Summary

Name	S22AA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q63ZE4
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAAIPGHRCW VHMLDNNTGS GNETGILSED ALLRISIPLD SNLRPEKCRR FVHPQWQLLH LNGTHSTSE ADTEPCVDGW VYDQSYFPST IVTKWDLVCD

image



Background

Alternative Names Solute carrier family 22 member 10 (Organic anion transporter 5)

Note

For Research Use Only , Not for Diagnostic Use