



Product Name: Recombinant Human P2RX1 Protein
Catalog #: PHH6069



Summary

Name	P2RX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P51575
Host	HEK293
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KQTQGYCAEH PEGGICKEDS GCTPGKAKRK AQGIRTGKCV AFNDTVKTCE IFGWCPVEVD DDIPRPALLR EAENFTLFIK NSISFPRFKV NRRNLVEEVN AAHMKTCLEFH KTLHPLCPVF QLGYYVQESG QNFSTLAEKG GVVGITIDWH CDLDWHVRHC RPIYEFHGLY EEKNLSPGFN FRFARHFVEN

image



Background

Alternative Names P2X purinoceptor 1 (P2X1) (ATP receptor) (Purinergic receptor)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human P2RX1 Protein
Catalog #: PHH6069

