


Summary

Name	SIDT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NXL6
Host	HEK293
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPGHPAKSPR QPPAPRRDPF DAARGADFDH VYSGVVNLST ENIYSFNYTS QPDQVTAVRV YVNSSSENLN YPVLVVVRQQ KEVLSWQVPL LFGQLYQRSY NYQEVSR TLC PSEATNETGP LQQLIFVDVA SMAPLGAQYK LLVTKLKH FQ LRTNVAFHFT ASPSQPYFL YKFPKD VDSV IIKVVSE MAY PCSVVS VQNI MCPVYDL DHN VEFNGVYQSM TKKAAITLQK KDFPGEQFFV VFVIKPEDYA CGGSFFIQEK ENQ TWNLQRK KNLEVTIVPS

image



Background

Alternative Names SID1 transmembrane family member 1

Product Name: Recombinant Human SIDT1 Protein
Catalog #: PHH6079



Note

For Research Use Only , Not for Diagnostic Use