



**Product Name: Recombinant Human SFRP5 Protein**  
**Catalog #: PHH6102**



## Summary

<b>Name</b>	SFRP5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-MBP & C-His
<b>Accession #</b>	Q5T4F7
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	57 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	IYPCRSLCEA VRAGCAPLME AYGFPWPEML HCHKFPLDND LCIAVQFGHL PATAPPVTKI CAQCMEHSA DGLMEQMCSS DFVVKMRIKE IKIENGDRKL IGAQKKKKLL KPGPLKRKDT KRLVLHMKNK

## image



## Background

**Alternative Names** Secreted frizzled-related protein 5 (sFRP-5) (Frizzled-related protein 1b) (FRP-1b) (Secreted apoptosis-related protein 3) (SARP-3)

## Note

For Research Use Only , Not for Diagnostic Use