



Product Name: Recombinant Human SPACA1 Protein
Catalog #: PHH6121



Summary

Name	SACA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9HBV2
Host	HEK293
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TENNDSETAE NYAPPETEDV SNRNVVKEVE FGMCTVTCGI GVREVILTNG CPGGESKCVV RVEECRGPTD CGWGKPISES LESVRLACIH

image



Background

Alternative Names Sperm acrosome membrane-associated protein 1 (Sperm acrosomal membrane-associated protein 32)

Note

For Research Use Only , Not for Diagnostic Use