


Product Name: Recombinant Human PROZ Protein
Catalog #: PHH6125



Summary

Name	PROZ_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P22891
Host	HEK293
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGNLEKECYE EICVYEEARE VFENEVVTDE FWRRYKGGSP CISQPCLHNG SCQDSIWGYT CTCSPGYEGS NCELAKNECH PERTDGCQHF CLPGQESYTC SCAQGYRLGE DHKQCVPHDQ CACGVLTSK RAPDLQDLPW QVKLTNSEGK

image



SDS



WB

Background

Alternative Names Vitamin K-dependent protein Z

Note

For Research Use Only , Not for Diagnostic Use