


Product Name: Recombinant Human SLC01A2 Protein
Catalog #: PHH6126



Summary

Name	SO1A2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P46721
Host	HEK293
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ENSSVVGINT SYEGIPQDLY VENDIFADCN VDCNCPSKIW DPVCGNNGLS YLSACLAGCE TSIGTGINMV FQNCSCIQTS GNSSAVLGLC DKGPDCSLML

image



Background

Alternative Names	Solute carrier organic anion transporter family member 1A2 (OATP-A) (Organic anion-transporting polypeptide 1) (OATP-1) (Sodium-independent organic anion transporter) (Solute carrier family 21 member 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use