


Product Name: Recombinant Human SLC01B7 Protein
Catalog #: PHH6128



Summary

Name	SO1B7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	G3V0H7
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LTLTYDGNP VRSHVDVPLS YCNSECNCDE SQWEPVCGNN GITYLSPCLA GCKSSSGNKE PIVFYNCSCV EVIGLQNKNY SAHLGECPRD DACTRKSYYV

image



Background

Alternative Names Putative solute carrier organic anion transporter family member 1B7 (Liver-specific organic anion transporter 3) (LST-3)

Note

For Research Use Only , Not for Diagnostic Use