



Product Name: Recombinant Human SLC01C1 Protein
Catalog #: PHH6129



Summary

Name	SO1C1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9NYB5
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DVAGLTVSYQ GTKPVSYHER ALFSDCNSRC KCSETKWEPM CGENGITYVS ACLAGCQTSN RSGKNIIFYN CTCVGIAASK SGNSSGIVGR CQKDNGCPQM

image



Background

Alternative Names	Solute carrier organic anion transporter family member 1C1 (Organic anion transporter F) (OATP-F) (Organic anion transporter polypeptide-related protein 5) (OAT-RP-5) (OATPRP5) (Organic anion-transporting polypeptide 14) (OATP-14) (Solute carrier family 21 member 14) (Thyroxine transporter)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SLC01C1 Protein
Catalog #: PHH6129

