


Product Name: Recombinant Human RNF215 Protein
Catalog #: PHH6141



Summary

Name	RN215_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9Y6U7
Host	HEK293
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GAGRGGARAV RVDVRLPRQD ALVLEGVRIG SEADPAPLLG GRLLLMDIVD AEQEAPVEGW IAVAYVGKEQ AAQFHQENKG SGPQAYPKAL VQQMRRALFL GASALLLIL NHNVVRELDI SQLLRPVIV LHYSSNVTKL LDALLQRTQA TAEITSGESL SANIEWKLTL WTCGLSKDG YGGWQDLVCL GGSRAQEQKP

image



Background

Alternative Names RING finger protein 215

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RNF215 Protein
Catalog #: PHH6141

